



SALSA m²

THE NEXT EVOLUTION
IN PERFORMANCE
AND STYLE





CHOOSE PERFORMANCE INDOORS...

Discover the next evolution in performance and style – the Quickie Salsa M².

Evolved from the Salsa M – the market leading mid-wheel-drive powerchair – the Salsa M² effortlessly balances superb indoor manoeuvrability with outstanding outdoor performance and head-turning style.

- ▶ Superior indoor manoeuvrability, thanks to the Salsa M²'s ultra-compact turning radius and Quickie's patented mid-wheel-drive technology.
- ▶ Get great access and control through narrow doorways, in lifts and around confined spaces with the ultra-compact base dimensions.
- ▶ Easily transfer into and out of your Salsa M², effortlessly access vehicles and under tables with the lower than ever before seat to floor height from 43cm (17") with powered COG tilt, or 45cm (17.5") with lift/tilt.



...AND OUT!

Evolved from highly-advanced technology and options from the Salsa M, the Salsa M² blends proven patented technology with a series of high-performance options from the Jive M.

- ▶ Enjoy a smooth and stable ride, even over rough terrain, with Quickie's patented all-wheel suspension - six independently suspended wheels absorb all those bumps.
- ▶ Quickie's anti-pitch technology eliminates pitching on slopes, and with all your weight located directly above the drive wheel to maximise traction you've never been in a safer driving position.
- ▶ More miles per charge! Take advantage of the ultra-efficient 4mph speed option powered by 60 Ah batteries and delivering 10% more charge than the original Salsa.
- ▶ Choose speed, torque and outright performance with the 120 Ah controller option and 6 mph, four-pole high torque motors.

Whatever your performance or practical requirements the Salsa M²'s range of options gives you the optimum balance between an exhilarating drive and great efficiency.

60 Ah Batteries

Range of 26 km (32 km with 60 Ah) according to the ISO standards.

2 or 4 Pole Motors

Performance motors with 6 and 10 kph options.

Anti-pitch & Suspension

Enjoy a smooth ride over uneven terrain with minimal body movement.



CHOOSE HEAD- TURNING STYLE...

The Salsa M² blends Quickie's true hallmarks of beautiful sleek design with an assertive stance and sporting performance - sure to turn heads.

- ▶ Stand out from the crowd with a selection of four body colours: matt black, classic white, vibrant red, or electric blue – there's a colour to suit every occasion and taste.
- ▶ Be seen with the Salsa M²'s LED bright lights and indicators integrated into the rear body panel for subtle styling.
- ▶ Make it your own! Customise your Salsa M² beyond the range of performance and style options with our Built-4-Me service, exclusively available through Quickie.

Whether you want to stand out from the crowd, turn heads in awe, or simply complement your outfit - with the Salsa M² you're always in style.

→ Four body colour options

built4me



...AND FUNCTION

Unparalleled performance and head-turning style doesn't mean compromising on practicality. The Salsa M² delivers great functionality:

- ▶ Flip-back armrests provide easy access in and out of your Salsa M² and allow for close access to desks.
- ▶ Effortlessly climb 7 cm (3") kerbs and thresholds.
- ▶ A 30 cm lift option gives you easy access to bars, cupboards and worktops.
- ▶ Drive your vehicle from the comfort of your Salsa M² or enjoy hassle-free transportation and peace of mind with the approved DAHL docking solution. The Salsa M² is crash tested and passed to ISO 7176-19.

Performance, style and practicality - without compromise. Choose from a wide range of options to perfectly tailor your Salsa M² to you and your lifestyle. Live without limits!



Powered tilt module (0° to 30°) with centre of gravity adjustment to suit correct users balance.



→ Fold-down backrest

The lift-off backrest ensures a compact size for easy transportation and storage.



SALSA M²

Approved for transport in a vehicle

The Salsa M² is transit tested in accordance with ISO 7176-19 and approved for transportation in a vehicle. The tests are conducted with a four-point-tie-down system and the optional DAHL docking system for an even easier fixation.

Every Salsa M² comes prepared for the DAHL docking station as standard - making fitting easier than ever!



Configure the Salsa M² to your individual needs – choose a wide range of options.
For further information on product options please see the Salsa M² order form.



// POWERED TILT MODULE WITH CENTRE OF GRAVITY ADJUSTMENT

Offers tilt from 0° up to 30° with an additional weight shift mechanism to suit correct users balance.



// THREE DIFFERENT MOTORS

- 6 kph 2-pole
- 10 kph 2-pole
- 10 kph 4-pole high torque performance motor



// POWERED RECLINE BACK

Offers a range from -5° to 30°. Easily retrofittable if developing needs requires higher support.



// EASY TO SERVICE

Easy battery access without tools. Top shroud is fixed with Velcro.



// TRUNK SUPPORT LATERALS

Multi-adjustable laterals offer advanced positioning of the trunk.



// FULL RANGE OF CONTROL DEVICES AVAILABLE

- VR2 standard control (up to two actuators)
- Up to four actuator control through a separate switch box
- Dual control
- R-net optional controller
- Full range of special controls

TECHNICAL DATA

SALSA M² // Mid-wheel powered wheelchair

Max. user weight:	140 kg	Battery size:	60 Ah
Front & rear castor wheels:	7" solid	Max. dynamic stability:	8° (14%)
Drive wheels:	13" solid and pneumatic	Overall width:	61 cm (2-pole motors)
Seat width:	41 – 51 cm	Turning radius:	600mm (fixed centre mount hanger)
Seat depth:	41 – 51 cm	Max. kerb climbing:	7 cm
Max. speed:	6 kph and 10 kph	Colour choices:	4 colours - red, blue, white & matt black
Max. range (ISO 7176-4):	26 km (32 km with 60 Ah)	Transport (ISO 7176-19):	4 point tie-down & Dahl docking system
Battery box to ground clearance:	65mm	Chair weight:	125 to 176 kg



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news from Quickie through our Facebook page and YouTube channel:
www.facebook.com/QuickieWheelchairsUK
www.youtube.com/quickiewheelchairseu

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



Sunrise Medical Limited
 Thorns Road
 Brierley Hill
 West Midlands
 DY5 2LD

Tel: +44 (0) 845 605 66 88
 Fax: +44 (0) 845 605 66 89
www.SunriseMedical.co.uk